

# PEN2 Peptide

**Cat. No. GTX85064-PEP**

|                     |                          |
|---------------------|--------------------------|
| <b>Applications</b> | Neutralizing /Inhibition |
| <b>Species</b>      | Human                    |

**Package**  
50 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
|--------------------|----------------------|

|                          |                 |
|--------------------------|-----------------|
| Neutralizing /Inhibition | Assay dependent |
|--------------------------|-----------------|

Not tested in other applications.

## Properties

|                      |                                                                             |
|----------------------|-----------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                      |
| <b>Buffer</b>        | PBS, 0.1% BSA                                                               |
| <b>Preservative</b>  | 0.02% Sodium azide                                                          |
| <b>Storage</b>       | Store the peptide at -20°C, stable for one year.                            |
| <b>Concentration</b> | 0.20 mg/ml (Please refer to the vial label for the specific concentration.) |
| <b>Conjugation</b>   | Unconjugated                                                                |

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).