

CTRP2 antibody

Cat. No. GTX85071

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
ICC/IF	10 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 30 kDa. ([Note](#))

Product Note These proteins are often highly modified post-translationally and migrate in SDS-PAGE at positions other than their predicted size.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	CTRP2 antibody was raised against recombinant human CTRP2.
Purification	Purified by DEAE chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

**GTX85071 ICC/IF Image**


ICC/IF analysis of Jurkat cells using GTX85071 CTRP2 antibody.

Working concentration : 20 µg/ml

**GTX85071 WB Image**

WB analysis of (A) 5 ng, (B) 25 ng, and (C) 50 ng CTRP2 recombinant protein using GTX85071 CTRP2 antibody.

Working concentration : 1 µg/ml

**GTX85071 ICC/IF Image**

ICC/IF analysis of Jurkat cells using GTX85071 CTRP2 antibody.

Working concentration : 10 µg/ml



For full product information, images and publications, please visit our [website](#).