

# DAD1 antibody

**Cat. No. GTX85083**

|                     |                   |
|---------------------|-------------------|
| <b>Host</b>         | Rabbit            |
| <b>Clonality</b>    | Polyclonal        |
| <b>Isotype</b>      | IgG               |
| <b>Applications</b> | WB, ICC/IF, ELISA |
| <b>Reactivity</b>   | Human, Mouse      |

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 0.5 - 2 µg/mL        |
| ICC/IF             | 10 µg/mL             |
| ELISA              | Assay dependent      |

Not tested in other applications.

**Calculated MW** 12 kDa. ( [Note](#) )

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS                                                                                                                                                                                                                        |
| <b>Preservative</b>  | 0.02% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| <b>Immunogen</b>     | DAD1 antibody was raised against a peptide corresponding to 14 amino acids near the C-terminus of human DAD1. The immunogen is located within the last 50 amino acids of DAD1.                                             |
| <b>Purification</b>  | Purified by antigen-affinity chromatography                                                                                                                                                                                |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTX85083 ICC/IF Image

ICC/IF analysis of HepG2 cells using GTX85083 DAD1 antibody.  
Working concentration : 10 µg/ml



### GTX85083 WB Image

WB analysis of HepG2 cell lysate using GTX85083 DAD1 antibody.  
Working concentration : (A) 0.5, (B) 1, and (C) 2 µg/ml



### GTX85083 ICC/IF Image

ICC/IF analysis of HepG2 cells using GTX85083 DAD1 antibody.  
Working concentration : 10 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).