

Her2 / ErbB2 antibody

Cat. No. GTX85113

Host	Chicken
Clonality	Polyclonal
Isotype	IgY
Applications	WB, IHC-P, IP, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	2.5 µg/mL
IP	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 138 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	erbB-2 antibody was raised against a 19 amino acid synthetic peptide near the carboxy terminus of human erbB-2. The immunogen is located within the last 50 amino acids of erbB-2.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX85113 IHC-P Image

IHC-P analysis of human liver tissue using GTX85113 Her2 / ErbB2 antibody.
Working concentration : 20 µg/ml



GTX85113 IHC-P Image

IHC-P analysis of human liver tissue using GTX85113 Her2 / ErbB2 antibody.
Working concentration : 2.5 µg/ml



GTX85113 WB Image

WB analysis of rat liver tissue lysate using GTX85113 Her2 / ErbB2 antibody.
Working concentration : (A) 1 and (B) 2 µg/ml



For full product information, images and publications, please visit our [website](#).