

NOD4 antibody

Cat. No. GTX85160

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, ELISA
Reactivity	Human, Mouse

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
ICC/IF	10 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 205 kDa. ([Note](#))

Product Note NOD4 antibody is predicted to not cross-react with other members of the NOD family of proteins.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	NOD4 antibody was raised against an 18 amino acid synthetic peptide near the amino terminus of human NOD4. The immunogen is located within amino acids 20 - 70 of NOD4.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX85160 ICC/IF Image


ICC/IF analysis of HeLa cells using GTX85160 NOD4 antibody.

Working concentration : 10 µg/ml


GTX85160 WB Image

WB analysis of EL4 cell lysate using GTX85160 NOD4 antibody.

Working concentration : 1 µg/ml


GTX85160 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX85160 NOD4 antibody.

Working concentration : 20 µg/ml



For full product information, images and publications, please visit our [website](#).