

SCRN3 antibody

Cat. No. GTX85171

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, ELISA
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
ICC/IF	20 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 49 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	SCRN3 antibody was raised against a 15 amino acid synthetic peptide near the center of human SCRIN3. The immunogen is located within amino acids 230 - 280 of SCRIN3.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx85171 WB Image

WB analysis of human breast tissue lysate using GTx85171 SCRn3 antibody.


Working concentration : (A) 1 and (B) 2 µg/ml



GTx85171 ICC/IF Image

ICC/IF analysis of 293 cells using GTx85171 SCRn3 antibody.

Working concentration : 20 µg/ml



GTx85171 ICC/IF Image

ICC/IF analysis of 293 cells using GTx85171 SCRn3 antibody.

Working concentration : 20 µg/ml



For full product information, images and publications, please visit our [website](#).