

CCDC69 antibody

Cat. No. GTX85190

Host	Chicken
Clonality	Polyclonal
Isotype	IgY
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-2 µg/mL
IHC-P	5-20µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 35 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	CCDC69 antibody was raised against a 16 amino acid synthetic peptide near the amino terminus of human CCDC69.The immunogen is located within amino acids 150 - 200 of CCDC69.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX85190 IHC-P Image

IHC-P analysis of human lung tissue using GTX85190 CCDC69 antibody.
Working concentration : 5 µg/ml




GTX85190 IHC-P Image

IHC-P analysis of human lung tissue using GTX85190 CCDC69 antibody.
Working concentration : 20 µg/ml



GTX85190 WB Image

WB analysis of rat pancreas tissue lysate using GTX85190 CCDC69 antibody.
Loading : 15µg
Dilution : 1µg/ml



GTX85190 WB Image

WB analysis of mouse pancreas tissue lysate using GTX85190 CCDC69 antibody.
Loading : 15µg
Dilution : 1µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).