

ZC3H12C Peptide

Cat. No. GTX85195-PEP

Applications	Neutralizing/Inhibition
Species	Human

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

Neutralizing/Inhibition	Assay dependent
-------------------------	-----------------

Not tested in other applications.

Properties

Form Liquid

Buffer PBS, 0.1% BSA

Preservative 0.02% Sodium azide

Storage Store the peptide at -20°C, stable for one year.

Concentration 0.20 mg/ml (Please refer to the vial label for the specific concentration.)

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).