

# OCIAD2 antibody

**Cat. No. GTX85233**

|                     |                   |
|---------------------|-------------------|
| <b>Host</b>         | Rabbit            |
| <b>Clonality</b>    | Polyclonal        |
| <b>Isotype</b>      | IgG               |
| <b>Applications</b> | WB, ICC/IF, ELISA |
| <b>Reactivity</b>   | Human             |

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 0.5 - 1 µg/mL        |
| ICC/IF             | 5 µg/mL              |
| ELISA              | Assay dependent      |

Not tested in other applications.

|                      |                                                 |
|----------------------|-------------------------------------------------|
| <b>Calculated MW</b> | 17 kDa. ( <a href="#">Note</a> )                |
| <b>Product Note</b>  | This antibody is specific for OCIAD2 C-Terminus |

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS                                                                                                                                                                                                                        |
| <b>Preservative</b>  | 0.02% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| <b>Immunogen</b>     | OCIAD2 antibody was raised against an 18 amino acid synthetic peptide from near the carboxy terminus of human OCIAD2. The immunogen is located within the last 50 amino acids of OCIAD2.                                   |
| <b>Purification</b>  | Purified by antigen-affinity chromatography                                                                                                                                                                                |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


## DATA IMAGES



### GTX85233 ICC/IF Image

ICC/IF analysis of A549 cells using GTX85233 OCIAD2 antibody.

Working concentration : 5 µg/ml



### GTX85233 WB Image

WB analysis of A549 cell lysate using GTX85233 OCIAD2 antibody.

Working concentration : (A) 0.5 and (B) 1 µg/ml



For full product information, images and publications, please visit our [website](#).