

SYPL2 antibody

Cat. No. GTX85272

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
ICC/IF	Assay dependent
IHC-P	2.5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 30 kDa. ([Note](#))

Product Note At least two isoforms of SYPL2 are known to exist. SYPL2 antibody will not cross-react with SYPL1.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	SYPL2 antibody was raised against a 15 amino acid synthetic peptide near the carboxy terminus of human SYPL. The immunogen is located within the last 50 amino acids of SYPL2.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX85272 ICC/IF Image

ICC/IF analysis of human spleen cells using GTX85272 SYPL2 antibody.
Working concentration : 20 µg/ml



GTX85272 WB Image

WB analysis of human spleen tissue lysate in (A) the absence and (B) the presence of blocking peptide using GTX85272 SYPL2 antibody.
Working concentration : 1 µg/ml



GTX85272 IHC-P Image

IHC-P analysis of human spleen tissue using GTX85272 SYPL2 antibody.
Working concentration : 2.5 µg/ml



For full product information, images and publications, please visit our [website](#).