

TTBK1 antibody

Cat. No. GTX85349

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
IHC-P	2.5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 143 kDa. ([Note](#))

Product Note Multiple isoforms of TTBK1 are known to exist. This antibody is predicted to not cross-react with TTBK2.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	TTBK1 antibody was raised against a 16 amino acid synthetic peptide near the carboxy terminus of human TTBK1. The immunogen is located within amino acids 1190 - 1240 of TTBK1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX85349 IHC-P Image

IHC-P analysis of human brain tissue using GTX85349 TTBK1 antibody.
Working concentration : 20 µg/ml



GTX85349 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX85349 TTBK1 antibody.
Working concentration : 2.5 µg/ml



GTX85349 WB Image

WB analysis of Jurkat cell lysate using GTX85349 TTBK1 antibody.
Working concentration : 1 µg/ml



For full product information, images and publications, please visit our [website](#).