

# OCC1 antibody

**Cat. No. GTX85350**

|                     |                  |
|---------------------|------------------|
| <b>Host</b>         | Rabbit           |
| <b>Clonality</b>    | Polyclonal       |
| <b>Isotype</b>      | IgG              |
| <b>Applications</b> | WB, IHC-P, ELISA |
| <b>Reactivity</b>   | Human            |

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 2 µg/mL              |
| IHC-P              | 2.5 µg/mL            |
| ELISA              | Assay dependent      |

Not tested in other applications.

**Calculated MW** 6 kDa. ( [Note](#) )

**Product Note** This antibody is specific for OCC-1 N-Terminus

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS                                                                                                                                                                                                                        |
| <b>Preservative</b>  | 0.02% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| <b>Immunogen</b>     | OCC-1 antibody was raised against a 16 amino acid synthetic peptide near the amino terminus of human OCC-1. The immunogen is located within amino acids 10 - 60 of OCC-1.                                                  |
| <b>Purification</b>  | Purified by antigen-affinity chromatography                                                                                                                                                                                |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

## Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


## DATA IMAGES



### GTX85350 WB Image

WB analysis of human spleen tissue lysate using GTX85350 OCC1 antibody.


Working concentration : 2 µg/ml



### GTX85350 IHC-P Image

IHC-P analysis of human spleen tissue using GTX85350 OCC1 antibody.

Working concentration : 20 µg/ml



### GTX85350 IHC-P Image

IHC-P analysis of human spleen tissue using GTX85350 OCC1 antibody.

Working concentration : 2.5 µg/ml



For full product information, images and publications, please visit our [website](#).