

METTL7A antibody

Cat. No. GTX85371

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	2 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 28 kDa. ([Note](#))**Product Note** At least two isoforms of Mettl7A are known to exist. This Mettl7A antibody is predicted to not cross-react with Mettl7B.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	MettL7A antibody was raised against a 12 amino acid synthetic peptide near the carboxy terminus of human Mettl7A. The immunogen is located within the last 50 amino acids of Mettl7A.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX85371 WB Image

WB analysis of A-20 cell lysate using GTX85371 METTL7A antibody.

Working concentration : 2 µg/ml



For full product information, images and publications, please visit our [website](#).