

Cathelicidin antibody

Cat. No. GTX85418

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 19 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Cathelicidin antibody was raised against a 17 amino acid synthetic peptide near the carboxy terminus of the human Cathelicidin. The immunogen is located within amino acids 50 - 100 of Cathelicidin.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX85418 WB Image


WB analysis of human spleen tissue lysate using GTX85418 Cathelicidin antibody.
Working concentration : 1 µg/ml


GTX85418 IHC-P Image

IHC-P analysis of human spleen tissue using GTX85418 Cathelicidin antibody.
Working concentration : 5 µg/ml


GTX85418 IHC-P Image

IHC-P analysis of human spleen tissue using GTX85418 Cathelicidin antibody.
Working concentration : 5 µg/ml


GTX85418 IHC-P Image

IHC-P analysis of human spleen tissue using GTX85418 Cathelicidin antibody.
Working concentration : 20 µg/ml



For full product information, images and publications, please visit our [website](#).