

## FEZ2 Peptide

Cat. No. GTX85422-PEP

<b>Applications</b>	Neutralizing /Inhibition
<b>Species</b>	Human

Package  
50 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Neutralizing /Inhibition	Assay dependent

Not tested in other applications.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS, 0.1% BSA
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store the peptide at -20°C, stable for one year.
<b>Concentration</b>	0.20 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Conjugation</b>	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).