

FNBP1L antibody

Cat. No. GTX85425

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
IHC-P	2.5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	70 kDa. (Note)
Product Note	This antibody is specific for TOCA-1 C-Terminus

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	TOCA-1 antibody was raised against a 17 amino acid synthetic peptide near the carboxy terminus of the human TOCA-1. The immunogen is located within the last 50 amino acids of TOCA-1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx85425 WB Image


WB analysis of mouse brain tissue lysate in (A) the absence and (B) the presence of blocking peptide using GTx85425 FNBP1L antibody.

Working concentration : 1 µg/ml


GTx85425 IHC-P Image

IHC-P analysis of human brain tissue using GTx85425 FNBP1L antibody.

Working concentration : 20 µg/ml


GTx85425 IHC-P Image

IHC-P analysis of human brain tissue using GTx85425 FNBP1L antibody.

Working concentration : 2.5 µg/ml



For full product information, images and publications, please visit our [website](#).