

CCDC98 Peptide

Cat. No. GTX85433-PEP

Applications	Neutralizing /Inhibition
---------------------	--------------------------

Species	Human
----------------	-------

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

Recommended dilution

Neutralizing /Inhibition

Assay dependent

Not tested in other applications.

Properties

Form	Liquid
-------------	--------

Buffer	PBS, 0.1% BSA
---------------	---------------

Preservative	0.02% Sodium azide
---------------------	--------------------

Storage	Store the peptide at -20°C, stable for one year.
----------------	--

Concentration	0.20 mg/ml (Please refer to the vial label for the specific concentration.)
----------------------	---

Conjugation	Unconjugated
--------------------	--------------

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).