

Nogo A antibody

Cat. No. GTX85448

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5 - 1 µg/mL
IHC-P	2.5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	130 kDa. (Note)
Product Note	At least five isoforms of Nogo are known to exist; this antibody is specific for NogoA.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	NogoA antibody was raised against a 23 amino acid synthetic peptide from near the amino terminus of human NogoA. The immunogen is located within the first 50 amino acids of NogoA.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX85448 WB Image

WB analysis of mouse brain tissue lysate using GTX85448 Nogo A antibody.
Working concentration : (A) 0.5 and (B) 1 μ g/ml



GTX85448 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX85448 Nogo A antibody.
Working concentration : 20 μ g/ml



GTX85448 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX85448 Nogo A antibody.
Working concentration : 2.5 μ g/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).