

RPA Interacting Protein antibody

Cat. No. GTX85485

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5 - 1 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 25 kDa. ([Note](#))

Product Note This antibody is specific for RPA Interacting Protein N-Terminus

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	RPA Interacting Protein antibody was raised against a 17 amino acid synthetic peptide from near the carboxy terminus of human RPA Interacting Protein. The immunogen is located within the first 50 amino acids of RPA Interacting Protein.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX85485 WB Image

WB analysis of Jurkat cell lysate using GTX85485 RPA Interacting Protein antibody.

Working concentration : (A) 0.5 and (B) 1 μ g/ml



For full product information, images and publications, please visit our [website](#).