

DDX58 antibody

Cat. No. GTX85488

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

References (4)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5 - 2 µg/mL
ICC/IF	Assay dependent
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 107 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	DDX58 antibody was raised against human GST-tagged DDX58 protein.
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx85488 ICC/IF Image

ICC/IF analysis of MCF7 cells using GTx85488 DDX58 antibody.
Dilution : 20 µg/ml




GTx85488 IHC-P Image

IHC-P analysis of human testis tissue using GTx85488 DDX58 antibody.
Dilution : 20 µg/ml



GTx85488 WB Image

WB analysis of C2C12 cell lysate using GTx85488 DDX58 antibody.
Working concentration : (A) 0.5, (B) 1, and (C) 2 µg/ml



GTx85488 IHC-P Image

IHC-P analysis of human heart tissue using GTx85488 DDX58 antibody.
Working concentration : 5 µg/ml



For full product information, images and publications, please visit our [website](#).