

ARMER antibody

Cat. No. GTX85516

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IHC-Fr, ELISA
Reactivity	Human, Mouse, Rat

References (4)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5 - 2 µg/mL
IHC-P	2 µg/mL
IHC-Fr	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	23 kDa. (Note)
Product Note	This antibody is specific for ARMER C-Terminus

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	ARMER antibody was raised against a peptide corresponding to 15 amino acids near the C-terminus of human ARMER. The immunogen is located within the last 50 amino acids of ARMER.
Purification	Purified by ion exchange chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTx85516 IHC-P Image

IHC-P analysis of mouse intestine tissue using GTx85516 ARMER antibody.


Working concentration : 2 µg/ml



GTx85516 WB Image

WB analysis of mouse small intestine tissue lysate using GTx85516 ARMER antibody.

Working concentration : (A) 0.5, (B) 1, and (C) 2 µg/ml



GTx85516 IHC-P Image

IHC-P analysis of human bladder tissue using GTx85516 ARMER antibody.

Working concentration : 2 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).