

PHAP3 antibody

Cat. No. GTX85687

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA, IHC
Reactivity	Human, Mouse, Rat

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5 - 2 µg/ml
IHC-P	2.5 µg/ml
ELISA	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Calculated MW 31 kDa. ([Note](#))

Product Note PHAP III has no cross-reaction to PHAP I and PHAP I2a.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to amino acids close to C-terminus of human PHAP III. The immunogen is located within amino acids 200 - 250 of PHAP III.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX85687 IHC-P Image

IHC-P analysis of human brain, cortex tissue using GTX85687 PHAP3 antibody.

Antigen retrieval : Heat-induced antigen retrieval

Dilution : 5 µg/ml



GTX85687 IHC Image

IHC analysis of human small intestine tissue using GTX85687 PHAP3 antibody.

Dilution : 10µg/ml



For full product information, images and publications, please visit our [website](#).