

STEAP4 antibody

Cat. No. GTX85695

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	2 µg/ml
IHC-P	5 µg/ml
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	52 kDa. (Note)
Product Note	This STEAP4 antibody does not cross-react with other STEAP proteins.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A 14 amino acid synthetic peptide near the center of human STEAP4. The immunogen is located within amino acids 210 - 260 of STEAP4.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTx85695 WB Image

WB analysis of rat liver tissue lysate using GTx85695 STEAP4 antibody.

Dilution : 0.5 µg/mL (Lane A) and 1 µg/mL (Lane B)



GTx85695 IHC-P Image

IHC-P analysis of human heart tissue using GTx85695 STEAP4 antibody.

Antigen retrieval : Heat-induced antigen retrieval

Dilution : 5 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).