

RIF1 antibody

Cat. No. GTX85728

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse

References (2)

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:3000
IHC-P	5 µg/ml
ELISA	1:5000 - 1:40000

Note : Heat induced antigen retrieval in pH 6.0 citrate buffer

Not tested in other applications.

Calculated MW 266 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	20mM Potassium Phosphate, 150mM NaCl
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide from mouse RIF1.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx85728 WB Image

WB analysis of 25 ug of MEF whole cell lysate using GTx85728 RIF1 antibody.
Dilution : 1 µg/ml



GTx85728 IHC-P Image

IHC-P analysis of human kidney tissue using GTx85728 RIF1 antibody.
Antigen retrieval : Heat-induced antigen retrieval
Dilution : 5 µg/ml



GTx85728 IHC-P Image

IHC-P analysis of human testis tissue using GTx85728 RIF1 antibody.
Antigen retrieval : Heat-induced antigen retrieval
Dilution : 5 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).