

STX12 antibody [15G2]

Cat. No. GTX85751

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat, Dog, Hamster

Package

25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/ml
IHC-P	20 µg/ml

Note : Heat induced antigen retrieval in pH 6.0 citrate buffer

Not tested in other applications.

Calculated MW	31 kDa. (Note)
Product Note	Additional bands of 27 and 37 kD may also be recognized , and are believed to correspond to degradation products of syntaxin 13.

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant cytoplasmic domain of rat syntaxin 13.
Purification	Protein G purified
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX85751 IHC-P Image

IHC-P analysis of human placenta tissue using GTX85751 STX12 antibody [15G2].

Antigen retrieval : Heat-induced antigen retrieval

Dilution : 20 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).