

CAND2 antibody

Cat. No. GTX85760

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IP, ELISA
Reactivity	Human

Package

25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:1000
IHC-P	1:500
IP	Assay dependent
ELISA	1:2000 - 1:10000

Note : Heat induced antigen retrieval in pH 6.0 citrate buffer

Not tested in other applications.

Calculated MW 135 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Serum
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	85 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Amino acids 1130-1143 of Human CAND2/TIP120B (C-terminal) coupled to KLH.
Purification	Purified by delipidation and defibrination
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX85760 IHC-P Image

IHC-P analysis of human kidney tissue using GTX85760 CAND2 antibody.

Antigen retrieval : Heat-induced antigen retrieval

Dilution : 1:500



For full product information, images and publications, please visit our [website](#).