

TMC6 antibody

Cat. No. GTX85762

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/ml
IHC-P	5 µg/ml
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 90 kDa. ([Note](#))

Product Note This EVER1 antibody does not cross-react with EVER2.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A 14 amino acid synthetic peptide near the N-terminus of human EVER1. The immunogen is located within the first 50 amino acids of EVER1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX85762 IHC-P Image

IHC-P analysis of human testis tissue using GTX85762 TMC6 antibody.

Antigen retrieval : Heat-induced antigen retrieval

Dilution : 5 µg/ml



GTX85762 WB Image

WB analysis of A-20 cell lysate using GTX85762 TMC6 antibody.

Dilution : 1 µg/mL (Lane A) and 2 µg/mL (Lane B)



For full product information, images and publications, please visit our [website](#).