

TMP21 antibody

Cat. No. GTX85770

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	3 µg/ml
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 25 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A 19 amino acid synthetic peptide near the C-terminus of human TMP21. The immunogen is located within amino acids 160 - 210 of TMP21.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX85770 WB Image


WB analysis of Raji cell lysate using GTX85770 TMP21 antibody.

Dilution : 0.5 μ g/mL (Lane A), 1 μ g/mL (Lane B), and 2 μ g/mL (Lane C)


GTX85770 IHC-P Image

IHC-P analysis of human brain tissue using GTX85770 TMP21 antibody.


Dilution : 20 μ g/ml


GTX85770 IHC-P Image

IHC-P analysis of human prostate tissue using GTX85770 TMP21 antibody.

Antigen retrieval : Heat-induced antigen retrieval

Dilution : 3 μ g/ml


GTX85770 IHC-P Image

IHC-P analysis of human breast tissue using GTX85770 TMP21 antibody.

Antigen retrieval : Heat-induced antigen retrieval

Dilution : 3 μ g/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).