

GLUT2 antibody [199017]

Cat. No. GTX85825

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ICC/IF, IHC-P, FCM
Reactivity	Human

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	10 µg/ml
IHC-P	10 µg/ml
FCM	Assay dependent

Note : Heat induced antigen retrieval in pH 6.0 citrate buffer

Not tested in other applications.

Product Note Based on flow cytometric tests on transfected cells, this antibody has no cross-reactivity with hGlut1 or hGlut3.

Properties

Form	Liquid
Buffer	PBS, 5% Trehalose
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human Glut2-transfected NS0 cells.
Purification	Protein G purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX85825 IHC-P Image

IHC-P analysis of human liver tissue using GTX85825 GLUT2 antibody [199017].

Antigen retrieval : Heat-induced antigen retrieval

Dilution : 10 µg/ml



For full product information, images and publications, please visit our [website](#).