

Noxa antibody

Cat. No. GTX85836

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/ml
IHC-P	10 µg/ml

Not tested in other applications.

Calculated MW 12 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A mixture of synthetic peptides containing amino acids 51-66 and 75-90 of mouse Noxa.
Purification	Protein G purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX85836 IHC-P Image

IHC-P analysis of human brain, cortex tissue using GTX85836 Noxa antibody.

Antigen retrieval : Heat-induced antigen retrieval

Dilution : 10 µg/ml



GTX85836 WB Image

WB analysis of A) human, B) mouse, and C) rat thymus tissue using GTX85836 Noxa antibody.

Dilution : 2 µg/ml



For full product information, images and publications, please visit our [website](#).