

LPA antibody

Cat. No. GTX85859

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P, ELISA
Reactivity	Human

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	20 µg/ml
ELISA	Assay dependent

Note : Heat induced antigen retrieval in pH 6.0 citrate buffer

Not tested in other applications.

Product Note No Cross-reactivity with Plasminogen or LDL.

Properties

Form	Liquid
Buffer	PBS, 0.2% BSA
Preservative	0.09% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human Lp(a).
Purification	Purified IgG
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX85859 IHC-P Image

IHC-P analysis of human liver tissue using GTX85859 LPA antibody.

Antigen retrieval : Heat-induced antigen retrieval

Dilution : 20 µg/ml



For full product information, images and publications, please visit our [website](#).