

CDH13 antibody

Cat. No. GTX85891

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, ELISA
Reactivity	Human, Mouse

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5-1 µg/mL
ICC/IF	20 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 78 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A 15 amino acid synthetic peptide near the N-terminus of human T-cadherin. The immunogen is located within amino acids 150 - 200 of T-cadherin.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTx85891 WB Image

WB analysis of 3T3 lysate using GTx85891 CDH13 antibody.
Dilution : 0.5 µg/mL (Lane A) and 1 µg/mL (Lane B)



GTx85891 ICC/IF Image

ICC/IF analysis of K562 cells using GTx85891 CDH13 antibody.
Dilution : 20µg/ml



GTx85891 IHC-P Image

IHC-P analysis of mouse brain tissue using GTx85891 CDH13 antibody.
Dilution : 5µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).