

Thyroid Hormone Receptor beta (cleaved Arg327) antibody

Cat. No. GTX86914

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500~1:1000

Not tested in other applications.

Calculated MW 53 kDa. ([Note](#))**Product Note** Thyroid Hormone Receptor beta (cleaved Arg327) antibody detects endogenous levels of fragment of activated THRB resulting from cleavage adjacent to Arg327.

Properties


Form	Liquid
Buffer	PBS, 150mM NaCl, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	The antiserum was produced against synthesized peptide derived from human Thyroid Hormone Receptor beta (278-327).
Purification	Purified by antigen-affinity chromatography From serum
Conjugation	Unconjugated
Note	For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 2

DATA IMAGES



GTx86914 WB Image

WB analysis of Jurkat cells treated with etoposide 25uM 24H using GTx86914 Thyroid Hormone Receptor beta (cleaved Arg327) antibody.



For full product information, images and publications, please visit our [website](#).