

## LPCAT2 antibody

Cat. No. GTX86989

|              |            |
|--------------|------------|
| Host         | Rabbit     |
| Clonality    | Polyclonal |
| Isotype      | IgG        |
| Applications | WB         |
| Reactivity   | Human      |

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution                | Recommended dilution |
|-----------------------------------|----------------------|
| WB                                | 1:500~1:1000         |
| Not tested in other applications. |                      |

Calculated MW 60 kDa. ( [Note](#) )

## Properties


|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 150mM NaCl, 50% Glycerol                                                                                                                                                                                              |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | The antiserum was produced against synthesized peptide derived from human LPCAT2 (321-370).                                                                                                                                |
| Purification  | Purified by antigen-affinity chromatography<br>From serum                                                                                                                                                                  |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |
| Note          | For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.                                                                                   |



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 2

## DATA IMAGES



### GTX86989 WB Image

WB analysis of K562 and HT-29 cell lysates using GTX86989 LPCAT2 antibody. The lane on the right is blocked with the synthesized peptide.



For full product information, images and publications, please visit our [website](#).