

SLC17A2 antibody

Cat. No. GTX87203

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2000
IHC-P	1:50-1:300

Not tested in other applications.

Calculated MW 47 kDa. ([Note](#))


Properties

Form	Liquid
Buffer	PBS, 150mM NaCl, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The antiserum was produced against synthesized peptide derived from human SLC17A2 (213-262).
Purification	Purified by antigen-affinity chromatography From serum
Conjugation	Unconjugated
Note	For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX87203 WB Image

WB analysis of HepG2 cell lysates using GTX87203 SLC17A2 antibody. The lane on the right is blocked with the synthesized peptide.




GTX87203 IHC-P Image

IHC-P analysis of human gastric adenocarcinoma tissue using GTX87203 SLC17A2 antibody
Antigen retrieval : Tris-EDTA, pH 9.0
Dilution : 1:200



GTX87203 WB Image

WB analysis of RAW264.7 whole cell lysate using GTX87203 SLC17A2 antibody.



GTX87203 WB Image

WB analysis of HepG2 whole cell lysate using GTX87203 SLC17A2 antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).