

# FFAR3 antibody

**Cat. No. GTX87383**

|                     |                |
|---------------------|----------------|
| <b>Host</b>         | Rabbit         |
| <b>Clonality</b>    | Polyclonal     |
| <b>Isotype</b>      | IgG            |
| <b>Applications</b> | ICC/IF, IHC-Fr |
| <b>Reactivity</b>   | Human, Mouse   |

References ( 1 )

Package

100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| ICC/IF             | 1:100~1:500          |
| IHC-Fr             | Assay dependent      |


Not tested in other applications.

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS, 150mM NaCl, 50% Glycerol                                                                                                                                                                                              |
| <b>Preservative</b>  | 0.02% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| <b>Immunogen</b>     | The antiserum was produced against synthesized peptide derived from human FFAR3 (11-60).                                                                                                                                   |
| <b>Purification</b>  | Purified by antigen-affinity chromatography<br>From serum                                                                                                                                                                  |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |
| <b>Note</b>          | For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.                                                                                   |



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX87383 ICC/IF Image**

ICC/IF analysis of LOVO cells using GTX87383 FFAR3 antibody. The picture on the right is blocked with the synthesized peptide.



For full product information, images and publications, please visit our [website](#).