

# ARHGAP22 antibody

**Cat. No. GTX87597**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Rabbit     |
| <b>Clonality</b>    | Polyclonal |
| <b>Isotype</b>      | IgG        |
| <b>Applications</b> | WB         |
| <b>Reactivity</b>   | Human      |

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500~1:1000         |

Not tested in other applications.

**Calculated MW** 77 kDa. ( [Note](#) )


## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS, 150mM NaCl, 50% Glycerol                                                                                                                                                                                              |
| <b>Preservative</b>  | 0.02% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| <b>Immunogen</b>     | The antiserum was produced against synthesized peptide derived from human ARHGAP22 (565-614).                                                                                                                              |
| <b>Purification</b>  | Purified by antigen-affinity chromatography<br>From serum                                                                                                                                                                  |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |
| <b>Note</b>          | For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.                                                                                   |



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX87597 WB Image

WB analysis of K562 cell lysates using GTX87597 ARHGAP22 antibody. The lane on the right is blocked with the synthesized peptide.



For full product information, images and publications, please visit our [website](https://www.genetex.com).