

## Arylsulfatase D antibody

Cat. No. GTX87652

|              |               |
|--------------|---------------|
| Host         | Rabbit        |
| Clonality    | Polyclonal    |
| Isotype      | IgG           |
| Applications | ICC/IF, IHC-P |
| Reactivity   | Human         |

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| ICC/IF             | 1:100~1:500          |
| IHC-P              | 1:50~1:100           |

Not tested in other applications.

## Properties


|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 150mM NaCl, 50% Glycerol                                                                                                                                                                                              |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | The antiserum was produced against synthesized peptide derived from human Arylsulfatase D (331-380).                                                                                                                       |
| Purification  | Purified by antigen-affinity chromatography<br>From serum                                                                                                                                                                  |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |
| Note          | For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.                                                                                   |



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 07 Page 1 of 2

## DATA IMAGES



### GTX87652 IHC-P Image

IHC-P analysis of human lung carcinoma tissue using GTX87652 Arylsulfatase D antibody. The picture on the right is blocked with the synthesized peptide.



### GTX87652 ICC/IF Image

ICC/IF analysis of MCF7 cells using GTX87652 Arylsulfatase D antibody. The picture on the right is blocked with the synthesized peptide.



For full product information, images and publications, please visit our [website](#).