

## hnRNP L antibody

**Cat. No. GTX87723**

<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	WB, ICC/IF, IHC-P
<b>Reactivity</b>	Human

**Package**  
100 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500~1:1000
ICC/IF	1:100~1:500
IHC-P	1:50~1:100

Not tested in other applications.


**Calculated MW** 64 kDa. ( [Note](#) )

## Properties


<b>Form</b>	Liquid
<b>Buffer</b>	PBS, 150mM NaCl, 50% Glycerol
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human hnRNP L (61-110).
<b>Purification</b>	Purified by antigen-affinity chromatography From serum
<b>Conjugation</b>	Unconjugated
<b>Note</b>	For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.




For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX87723 ICC/IF Image**

ICC/IF analysis of HeLa cells using GTX87723 hnRNP L antibody. The picture on the right is blocked with the synthesized peptide.


**GTX87723 IHC-P Image**

IHC-P analysis of human colon carcinoma tissue using GTX87723 hnRNP L antibody. The picture on the right is blocked with the synthesized peptide.


**GTX87723 WB Image**

WB analysis of HeLa, Jurkat, and HUVEC cell lysates using GTX87723 hnRNP L antibody. The lane on the right is blocked with the synthesized peptide.



For full product information, images and publications, please visit our [website](https://www.genetex.com).