

## Histone H2BK12ac (acetyl Lys12) antibody

Cat. No. GTX88011

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Application	WB, ICC/IF, IHC-P
Reactivity	Human, Monkey

Package  
100 µg

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500~1:1000
ICC/IF	1:100~1:500
IHC-P	1:50~1:100


Not tested in other applications.

## PROPERTIES


Form	Liquid
Buffer	PBS, 150mM NaCl, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	The antiserum was produced against synthesized peptide derived from human Histone H2BK12ac (8-59).
Purification	Purified by antigen-affinity chromatography. Non-acetylated peptide was removed by chromatography using corresponding non-acetylated peptide. From serum
Conjugation	Unconjugated
Note	For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.




For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX88011 IHC-P Image**

IHC-P analysis of human lung carcinoma tissue using GTX88011 Histone H2BK12ac (acetyl Lys12) antibody. The picture on the right is blocked with the synthesized peptide.


**GTX88011 WB Image**

WB analysis of COS7 cells treated with TSA 400nM (24h) lysate using GTX88011 Histone H2BK12ac (acetyl Lys12) antibody. The lane on the right is blocked with the synthesized peptide.


**GTX88011 ICC/IF Image**

ICC/IF analysis of HeLa cells using GTX88011 Histone H2BK12ac (acetyl Lys12) antibody. The picture on the right is blocked with the synthesized peptide.



For full product information, images and publications, please visit our [website](https://www.genetex.com).