

# Sodium/Potassium ATPase beta 1 blocking peptide

**Cat. No. GTX88074-PEP**

**Species** Human

**Package**  
100 µg

## Applications

**Product Note** The peptide was used in the production of GTX88074.

## Properties

**Form** Lyophilized powder

**Buffer** Reconstitute with 200µl distilled water to obtain a 0.5mg/ml peptide solution. Lyophilized from 5% (v/v) acetonitrile/H<sub>2</sub>O.

**Preservative** No preservatives

**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Region/Sequence** C-KTEISFRPNDPKSYE

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).