

DIXDC1 antibody, Internal

Cat. No. GTX88233

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

WB	0.5-2µg/ml
----	------------

Not tested in other applications.

Calculated MW	77 kDa. (Note)
----------------------	----------------------------------

Product Note	This antibody is expected to recognize both isoforms (NP_001033043.1; NP_219493.1).
---------------------	---

Properties

Form	Liquid
-------------	--------

Buffer	TBS, 0.5% BSA
---------------	---------------

Preservative	0.02% Sodium azide
---------------------	--------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
----------------	--

Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
----------------------	---

Immunogen	Peptide with sequence C-DVSYHQVDLERE, from the internal region of the protein sequence according to NP_001033043.1; NP_219493.1.
------------------	--

Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
---------------------	---

Conjugation	Unconjugated
--------------------	--------------

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

**GTX88233 WB Image**

WB analysis of human skeletal muscle lysate using GTX88233 DIXDC1 antibody, Internal.

Dilution : 0.5µg/ml

Loading : 35µg protein in RIPA buffer



For full product information, images and publications, please visit our [website](#).