

SA2 antibody, Internal

Cat. No. GTX88304

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.1-0.3µg/ml

Not tested in other applications.

Calculated MW 141 kDa. ([Note](#))

Product Note

This antibody is expected to recognize both reported isoforms (NP_001036214.1; NP_006594.3). Reported variants NP_001036215.1 and NP_001036214.1 represent identical protein: Reported variants NP_001036216.1 and NP_006594.3 represent identical protein.

Properties

Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence C-SRGSTVRSKSKSPST, from the internal region of the protein sequence according to NP_001036214.1; NP_006594.3.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX88304 WB Image

WB analysis of mouse spleen (A), mouse thymus (B) and rat spleen (C) lysates using GTX88304 SA2 antibody, Internal.

Dilution : 0.3µg/ml

Loading : 35µg protein in RIPA buffer


GTX88304 WB Image

WB analysis of K562 lysate using GTX88304 SA2 antibody, Internal.

Dilution : 0.03µg/ml

Loading : 35µg protein in RIPA buffer



For full product information, images and publications, please visit our [website](https://www.genetex.com).