

ARL4D blocking peptide

Cat. No. GTX88455-PEP

Species	Human
---------	-------

Package	100 µg
---------	--------

Applications

Product Note	The peptide was used in the production of GTX88455.
--------------	---

Properties

Form	Lyophilized powder
------	--------------------

Buffer	Reconstitute with 200µl distilled water to obtain a 0.5mg/ml peptide solution. Lyophilized from 5% (v/v) acetonitrile/H ₂ O.
--------	---

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Region/Sequence	C-KQDQPGALSAAEVEKR
-----------------	--------------------

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	---

	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
--	---



For full product information, images and publications, please visit our [website](#).