

MT2 blocking peptide

Cat. No. GTX88515-PEP

Species	Human
----------------	-------

Package 100 µg

Applications

Product Note	The peptide was used in the production of GTX88515.
---------------------	---

Properties

Form	Lyophilized powder
-------------	--------------------

Buffer	Reconstitute with 200µl distilled water to obtain a 0.5mg/ml peptide solution. Lyophilized from 5% (v/v) acetonitrile/H ₂ O.
---------------	---

Preservative	No preservatives
---------------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
----------------	--

Region/Sequence	C-QDASKGSHAEGQLQSP
------------------------	--------------------

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	---

	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
--	---



For full product information, images and publications, please visit our [website](#).