

## Tankyrase 1 antibody, N-term

Cat. No. GTX88517

|              |            |
|--------------|------------|
| Host         | Goat       |
| Clonality    | Polyclonal |
| Isotype      | IgG        |
| Applications | WB         |
| Reactivity   | Human      |

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 0.5-1.5µg/ml         |

Not tested in other applications.

|               |                                               |
|---------------|-----------------------------------------------|
| Calculated MW | 142 kDa. ( <a href="#">Note</a> )             |
| Product Note  | Immunogen shows no cross reactivity to TANK2. |

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | TBS, 0.5% BSA                                                                                                                                                                                                              |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 0.50 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                |
| Immunogen     | Fusion protein containing 120aa amino acids from N Terminal region of the protein sequence according to NP_003738.2.                                                                                                       |
| Purification  | Purified by ammonium sulphate precipitation followed by antigen affinity chromatography                                                                                                                                    |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX88517 WB Image**

WB analysis of HEK293 overexpressing Tankyrase1 (lane 2) and TANK2 (lane 4) (mock transfection in first lane) using GTX88517 Tankyrase 1 antibody, N-term. Lower panel shows the same lysates probed for alpha-Actin to show protein levels.

Dilution : 0.5µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).