

Septin 2 antibody, N-term

Cat. No. GTX88523

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.05-0.2µg/ml
IHC-P	2.5µg/ml

Not tested in other applications.

Calculated MW 41 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence SKQQPTQFINPET-C, from the N terminus of the protein sequence according to NP_001008491.1.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx88523 WB Image

WB analysis of human testis lysate using GTx88523 Septin 2 antibody, N-term.

Dilution : 0.05µg/ml


Loading : 35µg protein in RIPA buffer


GTx88523 IHC-P Image

IHC-P analysis of human cortex using GTx88523 Septin 2 antibody, N-term.

Antigen retrieval : citrate buffer pH 6


Dilution : 2.5µg/ml


GTx88523 IHC-P Image

IHC-P analysis of human tonsil using GTx88523 Septin 2 antibody, N-term.

Antigen retrieval : citrate buffer pH 6

Dilution : 2.5µg/ml


GTx88523 IHC-P Image

IHC-P analysis of human breast using GTx88523 Septin 2 antibody, N-term.

Antigen retrieval : citrate buffer pH 6

Dilution : 2.5µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).