

# Cytokeratin 13 antibody, C-term

**Cat. No. GTX88529**

<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	WB, IHC-P
<b>Reactivity</b>	Human

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.1-0.3µg/ml
IHC-P	2-4µg/ml

**Note : Human Tonsil shows strong layered, filamentous staining.**

Not tested in other applications.

<b>Calculated MW</b>	50 kDa. ( <a href="#">Note</a> )
<b>Product Note</b>	This antibody is expected to recognize the reported isoform a (NP_002265.2) only.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	TBS, 0.5% BSA
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Peptide with sequence C-TSNASGRRTSDVRRP, from the C Terminus of the protein sequence according to NP_705694.2.
<b>Purification</b>	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
<b>Conjugation</b>	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

## Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES




### GTx88529 WB Image

WB analysis of human lung lysate using GTx88529 Cytokeratin 13 antibody, C-term.

Dilution : 0.1µg/ml

Loading : 35µg protein in RIPA buffer



### GTx88529 IHC-P Image

IHC-P analysis of human tonsil using GTx88529 Cytokeratin 13 antibody, C-term.

Antigen retrieval : citrate buffer pH 6

Dilution : 2.5µg/ml



For full product information, images and publications, please visit our [website](#).