

Neurokinin 1 Receptor antibody, Internal

Cat. No. GTX89067

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

References (1) Package 100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution	
WB	1-3µg/ml	
Not tested in other applications.		
Calculated MW	46 kDa. (<u>Note</u>)	
Product Note	This antibody is expected to recognise both reported isoforms (NP_001049.1 and NP_056542.1).	

Properties	
Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence C-EIPGDSSDRYHEQ, from the internal region of the protein sequence according to NP_001049.1*; NP_056542.1*.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 29 Page 1 of 2



DATA IMAGES



GTX89067 WB Image

 $WB\ analysis\ of\ human\ testes\ lysate\ using\ GTX89067\ Neurokinin\ 1\ Receptor\ antibody,\ Internal.$

Dilution : $1\mu g/ml$

Loading : $35\mu g$ protein in RIPA buffer



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 29 Page 2 of 2