

MBD2 (isoform 1) antibody, Internal

Cat. No. GTX89182

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.3-1µg/ml
IHC-P	2-4µg/ml

Note : Human Colon shows nuclear staining of enterocytes.

Not tested in other applications.

Calculated MW	43 kDa. (Note)
Product Note	This antibody is expected to recognize isoform 1 (NP_003918.1) only.

Properties

Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence C-RNDPLNQNKGPDLN, from the internal region of the protein sequence according to NP_003918.1.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx89182 IHC-P Image


IHC-P analysis of human skin tissue using GTx89182 MBD2 (isoform 1) antibody, Internal.

Antigen retrieval : Steamed antigen retrieval with citrate buffer pH 6

Dilution : 2.5µg/ml

Antigen retrieval : Steamed antigen retrieval with citrate buffer pH 6

Dilution : 2.5µg/ml


GTx89182 IHC-P Image


IHC-P analysis of human adrenal gland tissue using GTx89182 MBD2 (isoform 1) antibody, Internal.

Antigen retrieval : Steamed antigen retrieval with citrate buffer pH 6

Dilution : 2.5µg/ml

Antigen retrieval : Steamed antigen retrieval with citrate buffer pH 6


Dilution : 2.5µg/ml


GTx89182 IHC-P Image

IHC-P analysis of human colon using GTx89182 MBD2 (isoform 1) antibody, Internal.

Antigen retrieval : citrate buffer pH 6

Dilution : 2.5µg/ml


GTx89182 WB Image

WB analysis of Jurkat cell lysate using GTx89182 MBD2 (isoform 1) antibody, Internal.

Dilution : 0.3µg/ml

Loading : 35µg protein in RIPA buffer



For full product information, images and publications, please visit our [website](https://www.genetex.com).